



HRF SWAB Experiment
Technical Interchange Meeting

Human Research Facility

Nancy Wilson

01/08/04

Surface, Water and Air Biocharacterization (SWAB) Experiment

Technical Interchange Meeting

January 8, 2004



HRF SWAB Experiment Technical Interchange Meeting	Human Research Facility	
	Nancy Wilson	01/08/04

Purpose of TIM

- Action Item from Phase 0/I safety review to return to discuss levels of containment for the SWAB water collection bags and SWAB sample tubes.



HRF SWAB Experiment Technical Interchange Meeting

Human Research Facility

Nancy Wilson

01/08/04

Experiment Overview: Inflight

- Inflight samples will be collected within 48 hours of any launch vehicle docking with ISS.
 - Collect at least three (3) water samples
 - 1 liter of SRV-K hot water
 - 1 liter of SRV-K cold water
 - 1 liter of SVO-ZV potable water
 - Collect any free floating condensate
 - Collection of two surface samples for each available ISS module currently on orbit.
 - Sample of 1 m³ of air for each of the available ISS modules currently on orbit.



HRF SWAB Experiment Technical Interchange Meeting

Human Research Facility

Nancy Wilson

01/08/04

Experiment Hardware List

- Air Sampling:
 - Air Sampling Device (ASD)
 - ASD Battery Packs
 - ASD Filter unit
- Surface Sampling:
 - HRF SWAB Tube (Manufactured by Charm Sciences)
- Water Sampling:
 - Potable Adapter Probe Assembly (SRV-K hot/cold)
 - SVO-ZV Adapter Probe Assembly
 - SWAB Temporary Water Bag(s) and SWAB Final Water Bag(s)
 - Contingency Syringe (condensate collection)



HRF SWAB Experiment Technical Interchange Meeting	Human Research Facility	
	Nancy Wilson	01/08/04

Surface Sampling Hardware

- Surface samples will be obtained with a cotton swab and stored in a buffer for analysis upon return to Earth.
- The swab tube (swab in vial w/buffer) will be used for surface sampling. The buffer is used to limit the rapid growth or proliferation of unwanted microbes as well as prolong the sample life.
- **Buffer solution has been determined to be toxicity level 1 for eye irritation (critical hazard).**



HRF SWAB Experiment Technical Interchange Meeting


Human Research Facility

Nancy Wilson

01/08/04

Original Surface Sampling Hardware

- Hardware selected prior to tox assessment, based on previous CHeCS hardware using identical buffer solution.



SWAB Surface Sampling Tube



HRF SWAB Experiment Technical Interchange Meeting


Human Research Facility

Nancy Wilson

01/08/04

Current Hardware Plans

- 2 levels of containment must be provided for 1ml of buffer solution.
- PocketSwab Plus
 - Provided with buffer solution in sealed portion of swab, sealed with foil seal
 - Launched in sealed ziplock
 - Removed for sample collection, swab pre-moistened with water
 - Buffer seal is punctured as swab is returned to body.
 - Replaced in ziplock following collection





HRF SWAB Experiment Technical Interchange Meeting

Human Research Facility

Nancy Wilson

01/08/04

Water Sampling Hardware

- The 1 liter water sample will be obtained with a specialized water port adapter which attaches the 1 liter sample bag to the ISS dispenser ports. The two water dispensing ports currently on the ISS require two different adapters.
- Hardware:
 - Potable Water Adapter Probe - SRV-K (hot and cold) dispensing ports
 - SVO-ZV Adapter Probe - SVO-ZV dispensing port.
 - 2 Water bag system – One to collect water and another that contains the fixative that stores the water for transport back to earth.
- In the event there is free floating condensate in any of the ISS modules, a contingency syringe will be used to collect the sample and transfer it to one of the spare 1 liter sample bags.
- Water bag contains fixative powder that has been classified as **Tox level 2 for eye irritant and inhalation.**



HRF SWAB Experiment Technical Interchange Meeting


Human Research Facility

Nancy Wilson

01/08/04

Water Sampling Hardware (Continued)

- SWAB Temporary Water Bag(s): This bag contains no fixative agent and will be connected to the water supply via the different adapter probes. A mechanically activated valve and cap are included on the inlet port.





HRF SWAB Experiment Technical Interchange Meeting


Human Research Facility

Nancy Wilson

01/08/04

Water Sampling Hardware (Concluded)

- SWAB Final Water Bag(s): This bag contains the fixative agent (tablets or water-soluble film bags) and will be connected to the Temporary Water Bag(s). The Final Bag is stowed inside a waterproof Bitran Bag during the transfer. A one-way valve and cap are included on the inlet port.





HRF SWAB Experiment Technical Interchange Meeting

Human Research Facility

Nancy Wilson

01/08/04

Implementation

- Tox level 2 powder would require 3 levels of containment.
- Tox level reduced to 0 if powder is in tablet form
- Final water bag is contained in a ziplock bag prior to sample collection. (film pouch provides 3rd level of containment)
- Water transfer is performed inside outer ziplock – approximately 5 minutes required for transfer
- Tox level is 1 after addition of water
- Water bag is 1 level of containment, outer ziplock bag is second level of containment.