

L O C K H E E D M A R T I N




2400 NASA Road One
Houston, Texas 77058

Human Research Facility (HRF) Rack 2
Refrigerated Centrifuge (RC)
Preliminary Fracture Control Summary Report

P/N SED46117400-301

CONTENTS

<u>Section</u>	<u>Page</u>
<u>CONTENTS</u>	i
1.0 <u>INTRODUCTION</u>	1-1
2.0 <u>FRACTURE CONTROL ASSESSMENT</u>	2-1
2.1 Crew-Induced Loads Analysis	2-1
2.2 Fail-Safe Analysis	2-1
2.3 Safe Return Configuration Analysis	2-1
2.4 Containment Analysis	2-2
2.5 Depressurization/Repressurization Analysis	2-2
2.6 Fatigue Analysis	2-2
2.7 Pressure System	2-3
2.8 Rotating Equipments	2-4
3.0 <u>SUMMARY</u>	3-1
APPENDIX A <u>CONTAINMENT ANALYSIS</u>	
APPENDIX B <u>FATIGUE ANALYSIS</u>	

Prepared By	Shaji Ommen	 LOCKHEED MARTIN SPACE OPERATIONS <i>Mechanical Systems Analysis Department</i>	Date	File Name
Checked By	D. Van		01 – 08 – 01	
Title Refrigerated Centrifuge (RC) Fracture Control Summary Report			Drawing No.	
			SEG46117400 – 301	
			Report No.	MSAD-01-0128

1.0 INTRODUCTION

The refrigerated centrifuge is intended to provide a system of separation of biological samples based on differing sample densities in a controlled temperature environment. RC will be a Commercial-Off-The-Shelf (COTS) unit, repackaged into a 12 Panel Unit (PU) drawer. It consists of two main components: (1) the refrigeration system and (2) the rotor assembly. The refrigeration method utilized by the centrifuge is vapor compression cycle. Vapor compression system consists of four components: compressor, condenser, evaporator, and expansion device.

During launch, landing, and on-orbit operation, the RC is rack mounted in a 12 PU active drawer. The refrigerated centrifuge will be used to separate biological samples such as blood and saliva. The front panel of the 12 PU drawer will open to expose the centrifuge rotor. The crewmember will select a rotor as identified in the experimental protocol procedures. The rotor is removed and replaced with the use of an Allen wrench that is provided as part of the centrifuge system. The samples will be loaded, and the door will be closed. The controls will be set for the appropriate time, temperature, rotor speed, and ramp up and down speeds.

Fracture control is implemented based on the Fracture Control Plan for Human Research Facility Payload and Racks, LS-71010. Structural analysis is performed based on the requirements in the Structural Verification Plan for Human Research Facility Payloads and Racks, LS-71012. RC design load factors are obtained from SSP Interface Requirements Document, SSP 57000, Revision D. Analysis results are documented in the HRF Rack 2 Refrigerated Centrifuge (RC) Stress Analysis Report, LMSEAT 33513, dated 10/2000.

Prepared By	Shaji Ommen	<div>LOCKHEED MARTIN</div> <div>SPACE OPERATIONS</div> <div>Mechanical Systems Analysis Department</div>	Date	01 – 08 – 01	File Name
Checked By	D. Van		Drawing No.	SEG46117400 – 301	
Title	Refrigerated Centrifuge (RC) Fracture Control Summary Report			Report No.	MSAD-01-0128





Figure 1 Refrigerated Centrifuge Assembly

Prepared By	Shaji Ommen	 LOCKHEED MARTIN SPACE OPERATIONS <i>Mechanical Systems Analysis Department</i>	Date	File Name
Checked By	D. Van		01 – 08 – 01	
Title Refrigerated Centrifuge (RC) Fracture Control Summary Report			Drawing No.	
			SEG46117400 – 301	
			Report No.	
			MSAD-01-0128	

2.0 **FRACTURE CONTROL ASSESSMENT**

Refrigerated Centrifuge panels are made of 7075-T7351 Aluminum. These panels are connected together by National Aerospace Standard (NAS) fasteners. All these panels and the fasteners are analyzed for launch and landing loads and have positive margins of safety. These panels and fasteners are classified as fail-safe by engineering judgment/analysis.

2.1 **Crew-Induced Loads Analysis**

The RC Front Door and the Stowage Drawer Front Panel are determined to be the weakest and the most vulnerable part for kick load. A 125 lb force is applied over a 4 inch x 4 inch square area. Analysis results show positive margins of safety. A crew-induced load of 50 lb is applied on the Front Handles. Analysis result shows positive margin of safety.


2.2 **Fail-Safe Analysis**

There are seven failure modes considered for the fail-safe analysis. These are (1) front panel to rack, (2) slide to side panel, (3) floating plate to bottom panel, (4) motor mounting ring to support assembly, (5) power controller module to bottom panel, (6) compressor to floating plate, and (7) condenser to floating plate. For all the cases the highly loaded fastener is assumed to be failed and the load normally carried by this fastener is now redistributed to the other fasteners. In all the cases analysis results show positive margins with a factor safety of 1.0.

2.3 **Safe Return Configuration Analysis**

According to the interpretation letter TA-93-037, "Structural Integrity Following Mechanism Failure", mechanism failure(s) which result in limit load redistribution will require structural verification on the redistributed loads if the failed condition is credible. When two-failure tolerance cannot be implemented, structural verification of the redistributed load path is required and the 1.4 factor of safety on limit loads must be maintained.

For the RC, the two credible failures are (1) floating plate to bottom panel fastener failure and (2) motor mounting ring to support assembly fastener failure. During on orbit operation the fasteners that connect the floating plate to the bottom panel and the ones that connect the motor mounting ring to support assembly will be removed to allow the vibration isolators to be in effect. In case of landing these fasteners should be reinstalled to overcome the flexibilities of the isolators. Analysis results show positive margins of safety with a factor of safety of 1.4

Prepared By	Shaji Ommen	LOCKHEED MARTIN  SPACE OPERATIONS <i>Mechanical Systems Analysis Department</i>	Date	File Name
Checked By	D. Van		01 – 08 – 01	
Title Refrigerated Centrifuge (RC) Fracture Control Summary Report			Drawing No.	
			SEG46117400 – 301	
			Report No.	MSAD-01-0128

2.4 Containment Analysis

For containment analysis, it must be shown that no part can attain sufficient kinetic energy to escape a container, which completely encompasses the part in order to prevent any hazard to the STS/ISS or crew (Ref. SSP 52005B, Section 6.2.1.1). The heaviest mass inside the RC drawer is the compressor, which is about 23.0 lb. In Section 2.2 it is shown that the floating plate, motor mounting ring, power controller module, compressor, and condenser connections are fail-safe. Therefore no containment analysis is necessary for these items.


The next heaviest item is the rotor stowed in the stowage drawer. The heaviest rotor has a mass of approximately 4.5 lb. A containment analysis is performed by assuming that this rotor is released from the stowage drawer and impacts the stowage drawer front panel. The result shows the stowage drawer front panels are adequate enough to contain the heaviest rotor. Detailed analysis is shown in Appendix A.

2.5 Depressurization/Repressurization Analysis

A depressurization/repressurization analysis is performed based on the requirements in SSP 57000, Revision D, Paragraph 3.1.1.2. The maximum delta pressure experienced during the depressurization and is negligible (0.00042 psi). In other words, the leak area is sufficient to prevent significant delta pressure build-up during depressurization or repressurization event. If this differential pressure is applied to the largest area and the smallest thickness panel, which is the top panel, which has an area of 401 in², the force imparted is only 0.168 lb. As a result, it can be stated that the effect of depressurization/repressurization is negligible.

2.6 Fatigue Analysis

A fatigue analysis of the RC is performed using NASA/FLAGRO crack growth program. Goddard load spectrum which is embedded within the program is used for the analysis. For the RC, the highly loaded element is the front panel flange, which is connected to the other panels is considered for the fatigue analysis. This section is 0.13 inch thick, 2.251 inch wide and has 0.23 inch diameter fastener holes. An initial corner crack size of 0.005 inch is used in this analysis. A tensile stress of 6.862 ksi, a bending stress of 23.359 ksi, and a bearing stress of 14.214 ksi on the front panel flange is used in this analysis. The result of the fatigue analysis is shown in Appendix B. The final fatigue life of the RC after applying a scatter factor, 4.0 is 16 missions.

Prepared By	Shaji Ommen	 LOCKHEED MARTIN SPACE OPERATIONS <i>Mechanical Systems Analysis Department</i>	Date	File Name
Checked By	D. Van		01 – 08 – 01	
Title Refrigerated Centrifuge (RC) Fracture Control Summary Report			Drawing No.	
			SEG46117400 – 301	
			Report No.	
			MSAD-01-0128	

The compressor takes low pressure, low temperature refrigerant gas and compresses it to high pressure, high temperature gas. Reciprocating pistons intake vapor at low pressure, and compress the vapor before sending it to the discharge line. From the compressor, the hot, high-pressure gas travels through the discharge line into the condenser. The condenser is the part of the system where the heat is rejected by, as the name implies, condensation. As the hot gas refrigerant cools, drops of liquid refrigerant form within the coil. Eventually, when the gas reaches the end of the condenser, it has condensed completely, that is, only liquid refrigerant is present.

The purpose of the expansion device in a vapor-compression refrigeration cycle is to control the flow of refrigerant to the evaporator. As the refrigerant leaves the condenser, it has cooled, and condensed to liquid, but is still under high pressure. In order for the liquid to absorb the necessary heat in the evaporator, its pressure must be reduced, which is accomplished within the expansion device. The refrigerated centrifuge uses a capillary tube to accomplish this. The evaporator is the component of the cycle, which actually absorbs the heat from the conditioned space. The evaporator is similar in construction to the condenser, but its function is opposite. As the fluid leaves the expansion device, it is a cool liquid. As it passes through the evaporator, it picks up heat from the room, and evaporates into a gaseous form. This evaporation is what enables the refrigerant to absorb the heat energy from the room.

As the refrigerant leaves the evaporator, it is returned to the cooled, low-pressure state, and is sent back to the compressor to begin the cycle again. Under normal circumstances the refrigerant will not wear out, it will be reused again and again, changing its physical form, but not its chemical composition. The refrigerated centrifuge uses R404a as its working fluid. R404a is environmentally friendly and has been given a toxicity level of 0 (reference memo #458). The pressure matrix for the RC is shown in Table 1. The Leak-Before-Burst (LBB) calculations for the fluid lines will be included in the final report.

2.8 Rotating Equipments

There are three fans in the Refrigerated Centrifuge assembly. There are two fans, P/N SEG46116060-702, which are mounted to the rear panel are HRF common fans. These fans measure 4.69 x 4.69 inches in size and spin at 5200 rpm at 28 VDC. If the fan speed is linearly extrapolated out to determine the worst-case speed at 32 VDC, then the speed will be 5943 rpm. The third fan, P/N W2E142-CC13-16 is connected to the condenser is a tubeaxial fan. This fan has a size of 5.91 x 6.77 inches and rotates at 2700 rpm at 115 VAC. For all the fans the size is less than 8 inches and maximum speed is less than 8000 rpm. Therefore, a containment analysis is not necessary.

The rotor motor and the rotor will be addressed in the final fracture control summary report.


Prepared By	Shaji Ommen	 LOCKHEED MARTIN SPACE OPERATIONS <i>Mechanical Systems Analysis Department</i>	Date	01 - 08 - 01	File Name
Checked By	D. Van		Drawing No. SEG46117400 - 301		
Title			Report No.		
Refrigerated Centrifuge (RC) Fracture Control Summary Report			MSAD-01-0128		

Table 1 Refrigerated Centrifuge Pressure Matrix

Pressure Matrix						
After two-fault failure, maximum system pressure is still controlled to be NO MORE THAN 350 psi						
Tubes, All Copper						
Size	Usage	MDP (psig)	(Burst) Pressure (psig)	Testing Organization	Reqd. FOS	Actl. FOS
2 mm	Capillary Tube	350	12000	ASME Refrigeration Piping Handbook	4.0	34.3
3/16" (5 mm)	Discharge	350	10000	Tested by LM	4.0	28.6
1/4" (6 mm)	Suction	150	6900	Tested by LM	4.0	46.0
5/16" (8 mm)	Evaporator	150	1200	Tested by LM	4.0	8.0
3/8" (10 mm)	Condenser	350	1700	ASME Refrigeration Piping Handbook	4.0	4.9
Fittings	P. Transducer	350	1500	Based on manufacturer data	4.0	4.3
Pressure Transducers	Control	350	1000	Based on manufacturer data	2.5	2.9
Filter/Drier						
Part No.	Manufacturer	MDP (psig)	Burst Pressure (psig)	Test Org.	Reqd. FOS	Actl. FOS
C-032-S	Sporlan	350	>2500	UL	2.5	7.1
Compressor Housing						
Part No.	Manufacturer	MDP (psig)	Burst Pressure (psig)	Test Org.	Reqd. FOS	Actl. FOS
AE820AT-900-J7	Tecumseh	150	>400	LM	2.5	2.7

Prepared By	Shaji Ommen	LOCKHEED MARTIN SPACE OPERATIONS <i>Mechanical Systems Analysis Department</i>	Date	01 – 08 – 01	File Name
Checked By	D. Van		Drawing No.	SEG46117400 – 301	
Title			Refrigerated Centrifuge (RC) Fracture Control Summary Report		
			Report No.	MSAD-01-0128	

3.0 SUMMARY

There is no fracture critical item in the Refrigerated Centrifuge. There is no glass, sealed container, or structural composite. All the panels and the fasteners connecting the panels together are classified as fail-safe. The RC stowage drawer will act as a structural containment for the stowed mass. The fracture control dispositions for all the items are shown in Table 2. RC meets all the fracture control requirements and is qualified for flight in EXPRESS Rack.

Table 2 Fracture Control Disposition Summary

ITEM OR PART	PART NUMBER	MATERIAL	NUMBER & TYPE OF CONNECTIONS	SAFETY FACTOR	CLASSIFICATION	NOTE
Top Panel	SDG46114565-001	7075-T7351 Aluminum	38 NAS1102-3	36.28	NFC (FS)	1
Bottom Panel	SDG46117378-301	7075-T7351 Aluminum	38 NAS1102-3	2.396	NFC (FS)	1
Right Side Panel	SDG46117367-301	7075-T7351 Aluminum	48 NAS1102-3	3.804	NFC (FS)	1
Left Side Panel	SDG46117367-301	7075-T7351 Aluminum	48 NAS1102-3	4.229	NFC (FS)	1
Front Panel	SDG46117366-301	7075-T7351 Aluminum	38 NAS1102-3	2.223	NFC (FS)	1
Rear Panel	SDG46117369-301	7075-T7351 Aluminum	38 NAS1102-3	2.809	NFC (FS)	1
Front Door	SEG46117371-301	7075-T7351 Aluminum	2 Latches & 2 Hinges	25.33	NFC (FS)	1
Stowage Drawer	SEG46117716-301	7075-T7351 Aluminum	63 NAS1101-08	5.532	NFC (FS)	1
Floating Plate	SDG46117388-301	7075-T7351 Aluminum	4 NAS1351N4	22.16	NFC (FS)	2
Support Assembly	SDG46117386-801	7075-T7351 Aluminum	8 NAS1101-3	5.194	NFC (FS)	1
Motor Mounting Ring	SDG46117391-301	7075-T7351 Aluminum	3 NAS1351N4	11.86	NFC (FS)	2
Slide	SEG46117459-301	7075-T7351 Aluminum	12 NAS1101-3	5.298	NFC (FS)	2
Internal Components					Contained	3
Pressure System					Pressurized Component	4
Fans	SEG46116060-702 W2E142-CC13-16		2 1		Contained	5


Prepared By	Shaji Ommen	 LOCKHEED MARTIN SPACE OPERATIONS <i>Mechanical Systems Analysis Department</i>	Date	File Name
Checked By	D. Van		01 – 08 – 01	
Title			Report No.	
Refrigerated Centrifuge (RC) Fracture Control Summary Report			MSAD-01-0128	

Table 2 Fracture Control Disposition Summary *Concluded*

ITEM OR PART	PART NUMBER	MATERIAL	NUMBER OF CONNECTIONS	SAFETY FACTOR	CLASSIFICATION	NOTE
Top Panel Fasteners	NAS1102-3	CRES A286	38	1.838	NFC (FS)	1
Bottom Panel Fasteners	NAS1102-3	CRES A286	38	1.725	NFC (FS)	1
Right Side Panel Fasteners	NAS1102-3	CRES A286	48	1.688	NFC (FS)	1
Left Side Panel Fasteners	NAS1102-3	CRES A286	48	1.600	NFC (FS)	1
Front Panel Fasteners	NAS1102-3	CRES A286	38	1.600	NFC (FS)	1
Rear Panel Fasteners	NAS1102-3	CRES A286	38	1.788	NFC (FS)	1
Stowage Drawer Fasteners	NAS1102-3	CRES A286	63	1.463	NFC (FS)	1
Support Assembly Fasteners	NAS1102-3	CRES A286	8	1.788	NFC (FS)	1
Front Panel to Rack Bolt	NAS1351N4	Heat-Res. Steel	6	1.584	NFC (FS)	2
Motor Ring to Support Assy. Bolt	NAS1351N4	Heat-Res. Steel	3	1.610	NFC (FS)	2
Floating Plate to Bottom Panel Bolt	NAS1351N4	Heat-Res. Steel	4	1.520	NFC (FS)	2
Slide to Side Panel Bolt	NAS1101-3	CRES A286	12	1.800	NFC (FS)	2
Power Controller to Bottom Panel Bolt	NAS1101-3	CRES A286	4	1.762	NFC (FS)	2
Compressor to Floating Plate Bolt	NAS1351N4	Heat-Res. Steel	4	1.615	NFC (FS)	2
Condenser to Floating Plate Bolt	NAS1101-3	CRES A286	4	1.954	NFC (FS)	2

NOTE 1: These items are connected with multiple fasteners and there is redundant load path. Fail-safe by engineering judgment.

NOTE 2: Fail-safe analysis and safe return configuration analysis are performed.

NOTE 3: See the containment analysis in Appendix A.

NOTE 4: See the pressure system description in Section 2.7.

NOTE 5: These fans are less than 8 inches in size and are less than 8000 rpm in speed.

APPENDIX A
CONTAINMENT ANALYSIS

A.0 CONTAINMENT ANALYSIS

For containment analysis, it must be shown that no part can attain sufficient kinetic energy to escape a container, which completely encompasses the part in order to prevent any hazard to the STS/ISS or crew (Ref. SSP 52005B, Section 6.2.1.1). The heaviest mass inside the RC drawer is the compressor, which is about 23.0 lb. In Section 7.0 it is shown that the floating plate, motor mounting ring, power controller module, compressor, and condenser connections are fail-safe. Therefore no containment analysis is necessary for these items. The next heaviest item is the rotor stowed in the stowage drawer. The heaviest rotor has a mass of approximately 4.5 lb.

A containment analysis is performed by assuming that this rotor is released from the stowage drawer and impacts the stowage drawer front panel. The front panel is 0.0625 inch thick. The worst-case net cg acceleration of 4.7g from MPLM Delta Design Cycle Couple Loads Analysis is used in the analysis, which is more conservative than the 1.5g suggested in SSP 52005B. It is assumed that the rotor with a nose radius of 0.5 inch travels 18 inches before it reaches the front panel. Structural containment is demonstrated using the following equation commonly referred to as the "Punch" equation.

$$T_R = \left[\frac{1}{2} * V^2 * \frac{W}{g} * \frac{1}{\pi d Y S_w} \right]^{1/2}$$

Where,

T_R = Maximum required wall thickness (inches) of the container to prevent escape of the component/part

W = Weight (pound-force) of the detached piece or part to be contained

g = Gravitational acceleration (inch per second squared)

V = Velocity (inches per second) that may be attained by that piece or part

d = Minimum profile diameter (inches) of piece or part that will produce a shear load on the container wall before escape by any particular piece or part resulting from a structural failure

$Y S_w$ = Yield strength (pound per square inch) of the container wall material

Velocity attained by the payload is computed by the following equation:

$$V = \sqrt{2aS_d}$$

Where,

$$a = \text{Steady state acceleration} = 4.7 \text{ g}$$

$$S_d = \text{Maximum travel distance of the rotor} = 18 \text{ in}$$

$$V = \sqrt{2 * 4.7 * 32.2 * 12 * 18} = 255.693 \text{ in/sec}$$

The front panel is made of 7075-T7351 Aluminum alloy and has yield strength of 57000 psi.
Calculating the required thickness for containment:

$$T_R = \left[\frac{1}{2} * 255.693^2 * \frac{4.5}{32.2 * 12} * \frac{1}{\pi * 1.0 * 57000} \right]^{[1/2]}$$

$$T_R = 0.0461 \text{ in}$$

Comparing the computed required thickness (0.0461 in) to the actual thickness of the front panel (0.0625 in), the actual thickness is higher. As a result, no penetration will occur.

APPENDIX B
FATIGUE ANALYSIS

B.0 FATIGUE ANALYSIS

```
FATIGUE CRACK GROWTH ANALYSIS
-----
DATE: 10/11/00    TIME: 14:43:08
(Computed: NASGRO Version 3.0.5, Apr 2000.)
U.S. customary units [inches, ksi, ksi sqrt(in)]

PROBLEM TITLE
-----
RC FATIGUE LIFE ANALYSIS

Crack Growth Model: Non Interaction
Equation/Table      : NASGRO Equation

GEOMETRY
-----
MODEL: CC02-Corner crack from hole in plate (2D).

Plate Thickness, t  =    0.1300
Plate Width, W      =    2.2510
Hole Diameter, D    =    0.2300
Edge to Hole Ctr, B =    0.3300
Poisson s ratio     =    0.33

FLAW SIZE: (User specified)

a  (init.) =  0.5000E-02
c  (init.) =  0.5000E-02
a/c (init.) =   1.000

MATERIAL
-----

MATL 1: 7075-T7351
      Plt & Sht; L-T

Material Properties:

:Matl:  UTS :  YS :  Kle :  Klc :  Ak :  Bk :  Thk :  Kc :  Keac :
: No.:      :      :      :      :      :      :      :      :      :
:-----:-----:-----:-----:-----:-----:-----:-----:
:  1 :  71.0:  62.0:  41.0:  29.0: 1.00: 1.00:  0.130:  56.4:      :

:Matl:----- Crack Growth Eqn Constants -----:
: No.:      C      :  n :  p :  q :  DKO : Cth+ : Cth- : Rcl:Alpha:Smax/:
:      :      :      :      :      :      :      :      :      :SIGo :
:-----:-----:-----:-----:-----:-----:-----:-----:
:  1 : 0.348D-07:2.529:0.50:1.00:  2.60:  2.00: 0.10:0.70: 1.90: 0.30:
```


RC FATIGUE LIFE ANALYSIS
MODEL: CC02

FATIGUE SCHEDULE BLOCK INPUT TABLE

GODDARD SPECTRUM

[Note: Stress = Input Value * Scale Factor]

Stress Scaling Factors for Block Case: 1

Scale Factor for Stress S0: 0.68620E-01
Scale Factor for Stress S1: 0.23359
Scale Factor for Stress S3: 0.14214

Schedule info. was input manually

Total No. of Blocks in Schedule = 1

Block Number and Case Correspondences
Block Number Block Case No.
From - To
1 - 1 1

Stresses: Tension, bending or pin

RC FATIGUE LIFE ANALYSIS
MODEL: CC02

FATIGUE SCHEDULE BLOCK INPUT TABLE

GODDARD SPECTRUM

SINGLE DISTINCT BLOCK

S	:	M:	NUMBER	:	S0	:	S1	:
T	:	A:	OF	:		:		:
E	:	T:	FATIGUE	:		:		:
P	:	L:	CYCLES	:	(t1)	:	(t2)	:

1:	:	1:	2.00	:	-100.00	:	100.00	:
2:	:	1:	4.00	:	-90.00	:	90.00	:
3:	:	1:	8.00	:	-80.00	:	80.00	:
4:	:	1:	15.00	:	-70.00	:	70.00	:
5:	:	1:	49.00	:	-60.00	:	60.00	:
6:	:	1:	81.00	:	-50.00	:	50.00	:
7:	:	1:	178.00	:	-40.00	:	40.00	:
8:	:	1:	641.00	:	-30.00	:	30.00	:
9:	:	1:	3120.00	:	-20.00	:	20.00	:
10:	:	1:	3405.00	:	-10.00	:	10.00	:
11:	:	1:	5019.00	:	-7.00	:	7.00	:
12:	:	1:	28853.00	:	-5.00	:	5.00	:
13:	:	1:	91655.00	:	-3.00	:	3.00	:

S	:	M:	NUMBER	:	S3	:	S	:
T	:	A:	OF	:		:		:
E	:	T:	FATIGUE	:		:		:
P	:	L:	CYCLES	:	(t1)	:	(t2)	:

1:	:	1:	2.00	:	-100.00	:	100.00	:
2:	:	1:	4.00	:	-90.00	:	90.00	:
3:	:	1:	8.00	:	-80.00	:	80.00	:
4:	:	1:	15.00	:	-70.00	:	70.00	:
5:	:	1:	49.00	:	-60.00	:	60.00	:
6:	:	1:	81.00	:	-50.00	:	50.00	:
7:	:	1:	178.00	:	-40.00	:	40.00	:
8:	:	1:	641.00	:	-30.00	:	30.00	:
9:	:	1:	3120.00	:	-20.00	:	20.00	:
10:	:	1:	3405.00	:	-10.00	:	10.00	:
11:	:	1:	5019.00	:	-7.00	:	7.00	:
12:	:	1:	28853.00	:	-5.00	:	5.00	:
13:	:	1:	91655.00	:	-3.00	:	3.00	:

Environmental Crack Growth Check for Sustained Stresses
(Kmax less than Keac): NOT SET

RC FATIGUE LIFE ANALYSIS
MODEL: CC02

FATIGUE SCHEDULE BLOCK STRESS TABLE

GODDARD SPECTRUM

S	:	M:	NUMBER	:	S0	:	S1	:
T	:	A:	OF	:		:		:
E	:	T:	FATIGUE	:	(ksi)	:	(ksi)	:
P	:	L:	CYCLES	:	(t1) : (t2)	:	(t1) : (t2)	:

1:	:	1:	2.00	:	-6.86:	6.86:	-23.36:	23.36:
2:	:	1:	4.00	:	-6.18:	6.18:	-21.02:	21.02:
3:	:	1:	8.00	:	-5.49:	5.49:	-18.69:	18.69:
4:	:	1:	15.00	:	-4.80:	4.80:	-16.35:	16.35:
5:	:	1:	49.00	:	-4.12:	4.12:	-14.02:	14.02:
6:	:	1:	81.00	:	-3.43:	3.43:	-11.68:	11.68:
7:	:	1:	178.00	:	-2.74:	2.74:	-9.34:	9.34:
8:	:	1:	641.00	:	-2.06:	2.06:	-7.01:	7.01:
9:	:	1:	3120.00	:	-1.37:	1.37:	-4.67:	4.67:
10:	:	1:	3405.00	:	-0.69:	0.69:	-2.34:	2.34:
11:	:	1:	5019.00	:	-0.48:	0.48:	-1.64:	1.64:
12:	:	1:	28853.00	:	-0.34:	0.34:	-1.17:	1.17:
13:	:	1:	91655.00	:	-0.21:	0.21:	-0.70:	0.70:

S	:	M:	NUMBER	:	S3	:	S	:
T	:	A:	OF	:		:		:
E	:	T:	FATIGUE	:	(ksi)	:	(ksi)	:
P	:	L:	CYCLES	:	(t1) : (t2)	:	(t1) : (t2)	:

1:	:	1:	2.00	:	-14.21:	14.21:	-14.21:	14.21:
2:	:	1:	4.00	:	-12.79:	12.79:	-12.79:	12.79:
3:	:	1:	8.00	:	-11.37:	11.37:	-11.37:	11.37:
4:	:	1:	15.00	:	-9.95:	9.95:	-9.95:	9.95:
5:	:	1:	49.00	:	-8.53:	8.53:	-8.53:	8.53:
6:	:	1:	81.00	:	-7.11:	7.11:	-7.11:	7.11:
7:	:	1:	178.00	:	-5.69:	5.69:	-5.69:	5.69:
8:	:	1:	641.00	:	-4.26:	4.26:	-4.26:	4.26:
9:	:	1:	3120.00	:	-2.84:	2.84:	-2.84:	2.84:
10:	:	1:	3405.00	:	-1.42:	1.42:	-1.42:	1.42:
11:	:	1:	5019.00	:	-0.99:	0.99:	-0.99:	0.99:
12:	:	1:	28853.00	:	-0.71:	0.71:	-0.71:	0.71:
13:	:	1:	91655.00	:	-0.43:	0.43:	-0.43:	0.43:

Environmental Crack Growth Check for Sustained Stresses
(Kmax less than Keac): NOT SET

RC FATIGUE LIFE ANALYSIS
MODEL: CC02

ANALYSIS RESULTS:

Schdl	Block Step	Cycles	Final Flaw Size a	Flaw Size c	K max a-tip	c-tip
10	1	1330300.	0.00948	0.00959	12.032	12.059
20	1	2660600.	0.01948	0.01973	14.516	14.622
30	1	3990900.	0.03402	0.03534	15.717	16.231
40	1	5321200.	0.05123	0.05617	16.410	17.750
50	1	6651500.	0.07063	0.08509	16.973	19.912
60	1	7981800.	0.09175	0.13024	17.445	23.504

FINAL RESULTS:

Crack outside geometric bounds:

c = 0.2154 B - D/2 = 0.2150

at Cycle No. 2302.69

of Load Step No. 9 Description: None

of Block No. 1

of Schedule No. 67

Crack Sizes: a = 0.114792 , c = 0.215358 , a/c = 0.5330

Total Cycles = 8783260.7

Execution time (hh:mm:ss): 00:00:01.0

Note: this is elapsed wall-clock time, not CPU time!